



Exceptional Training Makes all the Difference

Premature Anne™

Premature Anne[™] is a realistically proportioned 25-week preterm infant developed in collaboration with the American Academy of Pediatrics. Premature Anne is designed to facilitate training of healthcare professionals in proper care and resuscitation techniques for preterm infants.

Preterm infants are at greater risk and require specialized care by a highly skilled team of health professionals. Premature Anne is designed to meet the specific Neonatal Resuscitation Program® training requirements, defined by the American Academy of Pediatrics, to help ensure better patient outcomes for the tiniest of babies.

Premature Anne[™]

Premature Anne[™] Task Trainer

- · Anatomically accurate, realistic airway
- ET tube insertion
- Sellick Manoeuvre
- Positive Pressure Ventilation
- · Right mainstem intubation
- Suctioning
- OG/NG tube insertion
- Bilateral and unilateral chest rise and fall with mechanical ventilation
- Unilateral chest movement (right mainstem intubation) with mechanical ventilation
- Realistic Compressions
- Patent, cuttable umbilicus with venous and arterial access for bolus or infusion
- Simulated blood flashback upon cannulation of umbilical vein
- Peripheral IV access
- Auscultation of lung sounds during ventilation

Premature Anne[™] with SimPad[®] PLUS

- Easy to use, intuitive SimPad Plus user interface with 8 pre-programmed scenarios written by the American Academy of Pediatrics and supports the Neonatal Resuscitation Program[®]
- Anatomically accurate, realistic airway
- ET tube insertion
- Sellick Manoeuvre
- Positive Pressure Ventilation
- Right mainstem intubation
- Suctioning
- OG/NG tube insertion
- Bilateral and unilateral chest rise and fall with mechanical ventilation
- Cyanosis
- Unilateral chest movement with mechanical ventilation
- Realistic Compressions
- Patent, cuttable umbilicus with venous and arterial access for bolus or infusion
- Simulated blood flashback upon cannulation of umbilical vein
- Peripheral IV access
- Auscultation of lung sounds during ventilation
- Heart sounds
- Vocal sounds

ORDERING INFORMATION

290-00050 Premature Anne Task Trainer

Includes: Manikin, Hat/Blanket, Umbilical Cords (3), SimMom Birthing Lubricant, Airway Lubricant, 3M Tape, 42mm Facemask, 10ml Syringe, Umbilical Blood Concentrate IPI and User Guide.

290-00150 Premature Anne Task Trainer Twin Pack

Each Includes: Manikin, Hat/Blanket, Umbilical Cords (3), SimMom Birthing Lubricant, Airway Lubricant, 3M Tape, 42mm Facemask, 10ml Syringe, Umbilical Blood Concentrate IPI and User Guide.

295-55101 Premature Anne with SimPad Plus

Includes: Premature Anne Manikin, SimPad PLUS, NRP Scenarios (8 pre-loaded), Premature Anne Software for SimPad Plus (pre-loaded), Hat/Blanket, Umbilical Cords (3), SimMom Birthing Lubricant, Airway Lubricant, 3M Tape, 42mm Facemask, I Oml Syringe, Umbilical Blood Concentrate IPI and User Guide.



Operation

Premature Anne manikin can be operated by a SimPad® PLUS remote and can utilize a simulated patient monitor. Incorporate SimCapture for a complete debriefing solution.

SimPad PLUS







Optimize Your Investment

Laerdal offers comprehensive Educational and Technical Services designed to support healthcare professionals successfully implement and sustain simulation training programs. Our team of educators and technicians have years of industry related experience, with the ability to provide a wide range of options to help you maintain your simulation equipment and gain the critical skills necessary to optimize performance and learning outcomes.

