



GD/F30S Gynaecological Examination Model

Features:

1. Palpation of normal and pregnant uterus
2. Bi-manual pelvic and tri-manual pelvic examination
3. Vaginal examination
4. Visual recognition of normal and abnormal cervixes
5. Contraceptive instrument insertion and removal
6. Observation of uterus, ovary, fallopian tubes, and ligament teres
7. Normal and abnormal uterus and attachment model
8. Place and remove IUD by using the guide fork of contraceptive ring
9. Gravid uterus (Five-mouth fetus)
10. Extra-uterine pregnancy (ampullary pregnancy)
11. Salpingemphraxis



GD/FW2 Maternity Examination Model

- Specifically designed for professional training. The various functions of this model are as follow:-
- Palpation, embryo heart sound auscultation, pelvic measurement & mammary nursing training;
- Inflating uterus through ball to make model to be lifelike;
- Adjustable heart sound speed, volume;
- Heart sound can be heard not only from stethoscope but also from microphone of front board of model.