



GD/F30S Gynaecological Examination Model

Features:

- 1.Palpation of normal and pregnant uterus
- 2.Bi-manual pelvic and tri-manual pelvic examination
- 3.Vaginal examination
- 4.Visual recognition of normal and abnormal cervixes
- 5.Contraceptive instrument insertion and removal
- 6.Observation of uterus, ovary, fallopian tubes, and ligament teres
- 7.Normal and abnormal uterus and attachment model
- 8.Place and remove IUD by using the guide fork of contraceptive ring
- 9.Graavid uterus(Five-mouth fetus)
- 10.Extra-uterine pregnancy(ampullary pregnancy)
- 11.Salpingemphraxis



GD/FW2 Maternity Examination Model

- Specifically deigned for professional training. The various functions of this model are as follow:-
- Palpation, embryo heart sound auscultation, pelvic measurement & mammary nursing training;
- Inflating uterus through ball to make model to be lifelike;
- Adjustable heart sound speed, volume;
- Heart sound can be heard not only from stethoscope but also from microphone of front board of model.