



GD/F3B Vulva Suturing Training Simulator

Right Posterior incision ,left Posterior incision or median incision can be selected to exercise perineum cutting and suturing.



GD/F8 Vacuum Delivery Model

Features:

The simulator can help to demonstrate all the process of delivery. The life-size pelvic cavity has main anatomical structures.

It can help to train the skills of delivery. Smooth, soft skull with fontanelles for realistic vacuum delivery.

GD/F6 Leopold Maneuvers Simulator

Features:

1. Vertical section model, designed in accordance with the human physiological anatomy structure.
2. Accurate anatomic marks, obvious outline of uterus and bony pelvis.
3. Delivery mechanical parts demonstrate the whole delivery process, including engagement, descending, flexion, internal rotation, extension, resoration, external rotation and shoulder delivery.
4. Simulational birth canal skin, simulate the real delivery situation.
5. The fetal skin is soft and fontanel is visible.

