AIRWAY MANAGEMENT TRAINER



NG Tube and Trach Care Trainer

Head features anatomical landmarks, trachea, esophagus, simulated lungs, and stomach Lungs and stomach may be filled with fluid for realistic practice of many procedures:

- Tracheostomy care
- Tracheal suctioning
- •NG tube insertion and removal
- •NG tube irrigation, instillation, and monitoring
- •Feeding tube insertion and removal
- Gastric lavage and gavage
- •Nasoenteric and esophageal tube insertion, care, and removal
- Oropharyngeal and nasopharyngeal insertion and suctioning
- •Insertion, securing, and care of endotracheal tubes
- •375-10001 NG Tube and Trach Care Trainer



Deluxe Difficult Airway Trainer

The Deluxe Difficult Airway Trainer features a manually inflatable tongue to simulate obstructed airway and is designed for training the management of difficult airways.

- Manually inflatable tongue to simulate obstructed airway or tongue edema
- •Anatomically correct internal anatomy facilitates airway management skills using manual maneuvers and mechanical devices
- •Updated airway allows for insertion of standard airway devices with the addition of the LMA, Combitube and King LT
- Oropharyngeal and nasopharyngeal airway insertion
- Endotracheal and right mainstem intubation
- Needle or surgical cricothyrotomy
- •Retrograde intubation
- •Bag-Valve-Mask ventilation
- Tracheal suctioning
- •261-10001 Deluxe Difficult Airway Trainer